

# Canon

image*ANYWARE*

## Program Imaging Mode of Canon CXDI DR System Software

### Optimizing Digital Radiography Workflow

The Program Imaging Mode of Canon CXDI DR system software assists in optimizing Digital Radiography workflow. This high performance solution integrates the RIS with the Canon CXDI DR system software, utilizing the Procedure Code ID of the RIS to configure the exam. This feature increases the efficiency of the X-ray technologist by allowing improved workflow when preparing for the exam and the selection of views.



The flexibility of Canon portable DR Systems

### Other Advantages of the Program Imaging Mode Include:

- Helps to reduce potential errors in the selection of radiographic views
- Promotes efficiency of patient exams by displaying only required views
- Customize sequence of image acquisition
- Add, edit or delete views as required
- User-friendly: Ideal for beginning technologists, part-time technologists and/or students

#### RIS

Patient Name: Doe, John

Patient ID: 86088644

Description: Knee 4 Views

Accession Number: 34665687

Procedure Code ID

Worklist



Program Imaging Mode:  
Eliminates Unnecessary Buttons




Manual Imaging Mode:  
Populates With All Buttons




## Workflow Enhancement

Program Imaging Mode enhances workflow by eliminating many of the steps required to perform patient exams. Automatically assigns accession number to view in multi-accession exams.



*Performed views are grayed out and subsequent view is automatically selected.*



## Workflow Automation

The exposure preparation to obtain subsequent views on a patient is automatically selected by the CXDI system. In addition, performed views are grayed out so as to distinguish completed views.



## Workflow Flexibility

In trauma situations when the patient information is not obtained through the worklist, information and views can be added, as required.

*Add, edit or delete views as required.*



## CANON MEDICAL SYSTEMS

A Division of Canon U.S.A., Inc.

15955 Alton Parkway

Irvine, CA 92618-3731 U.S.A.

Tel: (800) 970-7227 • Fax: (949) 753-4164

<http://www.usa.canon.com/DR>